



# Recombinant *Drosophila melanogaster* 26S protease regulatory subunit 4 (Pros26.4)

<b>Product Code</b>	CSB-EP018888DLU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P48601
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGQNQSAQGG AGEKKDDKDK KKKYEPPIPT RVGKKKRRAK GPDAAMKLPQ VTPHTRCRLK LLKLERIKDY LMMEDEFIRN QERLKPQDEK NEEERSKVDD LRGTPMSVGN LEEIIDDNHA IVSTSVGSEH YVSILSFVDK DQLEPGCSVL LNHKVHAVVG VLSDDTDPMV TVMKLEKAPQ ETYADIGGLD TQIQEIKESV ELPLTHPEYY EEMGIKPPKG VILYGPPGTG KTLAKAVAN QTSATFLRVV GSELIQKYLK DGPKLVRELF RVAEEHAPSI VFIDEIDAVG TKRYDSNSGG EREIQRMTLE LLNQLDGFDS RGDVKVIMAT NRIETLDPAL IRPGRIDRKI EFPLPDEKTK RRIFTIHTSR MTLAEDVNLS ELIMAKDDLK GADIKAICTE AGLMALRERR MKVTNEDFKK SKESVLYRKK EGTPEGLYL
<b>Source</b>	<i>E.coli</i>
<b>Target Names</b>	Rpt2
<b>Protein Names</b>	Recommended name: 26S protease regulatory subunit 4 Short name= P26s4
<b>Expression Region</b>	1-439
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.