



Recombinant Rat Matrix metalloproteinase-9 (Mmp9)

Product Code	CSB-BP014679RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50282
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	FQT FDGDLKWHHH NITYWIQSYT EDLPRDVIDD SFARAFVWS AVTPLTFTRV YGLEADIVIQ FGVAEHGDGY PFDGKDGLLA HAFPPGPGIQ GDAHFDDDEL WSLGKGAVVP TYFGNANGAP CHFPTTFEGR SYLSCTTDGR NDGKPWCGTT ADYDTRKYG FCPSENLYTE HGNGDGKPCV PPFIFEGHSY SACTTKGRSD GYRWCATTAN YDQDKADGFC PTRADVTVTG GNSAGEMCVF PFVFLGKQYS TCTSEGRSDG RLWCATTSNF DADKKWGFPC DQGYSFLVA AHEFGHALGL DHSSVPEALM YPMYHYHEDS PLHEDDIKGI HHLYGRGSKP DPRPPATTA EPQPTAPPTM CSTAPPMAYP TGGPTVAPTG APSPGPTGPP TAGPSEAPTE SSTPDDNPCN VDFVDAIADI QGALHFFKDG RYWKFSNHGG NQLQGPFLIA RTWPAFPSKL NSAFEDPQPK KIFFFLWAQM WVYTGGQSVLG PRSLDKLGLG SEVTLVTGLL PRRGGKALLI SRERIWKFDL KSQKVPQSV TRLDNEFSGV PWNSHNVFQY QDKAYFCHDK YFWRVVFHNR VNQVDHVAYV TYDLLQCP
Source	Baculovirus
Target Names	Mmp9
Protein Names	Recommended name: Matrix metalloproteinase-9 Short name= MMP-9 EC= 3.4.24.35 Alternative name(s): 92 kDa gelatinase 92 kDa type IV collagenase Gelatinase B Short name= GELB
Expression Region	108-708
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP s are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. The enzyme encoded by this gene degrades type IV and V collagens. Studies in rhesus monkeys suggest that the enzyme is involved in IL-8-induced mobilization of hematopoietic progenitor cells from bone marrow,



and murine studies suggest a role in tumor-associated tissue remodeling.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.