



Recombinant 40S ribosomal protein SA (rps-0)

Product Code	CSB-YP020485CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P46769
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	SGGAAHSAL TEEDVMKLLA TQAHLGSTNL NFQMQQYVYK RRFDGPNIIN VKKTWEKLLL AARAIAAVEN PADVVVVSAR PYAQRALLKF AAHTGATAIF GRFSPGCLTN QIQKTFKEPR LLVISDPRID HQAVTEASYV GVPVISFVNT ESPLKLIDIG VPCNNKGERS IGLMWWMLAR EILILRGKIS RQTGFVLEGK EIMPDLYFYR DPTETEKEET GAHADVAEAQ EYQQPTDIDF TTQGGKVDDW AAETATWTAE TKTTEEWANA PTQSNW
Source	Yeast
Target Names	rps-0
Protein Names	Recommended name: 40S ribosomal protein SA
Expression Region	2-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.