



Recombinant Drosophila melanogaster Calcineurin subunit B type 1 (CanB)

Product Code	CSB-MP018575DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P48451
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MGNETSLPMD MCSNFDADEI RRLGKRFRKL DLDNSGALSI DEFMSLPELQ QNPLVQRVID IFDADNGEV DFKEFIQGVV QFSVRGDKLS KLRFAFRIYD MDNDGYISNG ELFQVLKMMV GNNLKDTQLQ QIVDKTICFA DKDEDGKISF DEFCSVVGNT DIHKKMVVDV
Source	Mammalian cell
Target Names	CanB
Protein Names	Recommended name: Calcineurin subunit B type 1 Alternative name(s): Protein phosphatase 2B regulatory subunit 1
Expression Region	1-170
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.