



Recombinant Drosophila melanogaster MAP kinase-activated protein kinase 2 (MAPk-Ak2)

Product Code	CSB-MP013473DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49071
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MLSLQNQRQP KTTPLTDDYV TSNTVLGYGI NGKVVQCTHR RTQQNYALKV LLDSERARRE VDLHWRVSGC RYIVNIIDVY ENTFKDRKCL LVMMECMEGG ELFQRIQDKA DGAFTEREAA QIMHEICAAV DYLSHRDIAH RDLKPENLLY TTTQP NATLK LTDFGFAKET FTSYTLQTPC YTPYYVAPEV LGPEKYDKSC DIWSLGVVMY IIMCGFPFVY SNHGLAISPG MKNRIRTGQY DFPDPEWTNV SQA AKDLIKG MLNVDPSKRL RIQDVISNKW IAQYNAV PQT PLCTGRMLKE AEETWPEVQE EMTRSLATMR VDYDQM QIKA LDKSNNPLLT KRRKKIEEME LYMANATRN
Source	Mammalian cell
Target Names	MAPk-Ak2
Protein Names	Recommended name: MAP kinase-activated protein kinase 2 Short name= MAPK-activated protein kinase 2 Short name= MAPKAP kinase 2 Short name= MAPKAPK-2 EC= 2.7.11.1
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.