



Recombinant *Saccharomyces cerevisiae* DNA polymerase delta subunit 3 (POL32)

Product Code	CSB-YP018305SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P47110
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDQKASYFIN EKLFTVEKPV LFTDLIHHLK IGPSMAKKLM FDYYKQTTNA KYNCVVICCY KDQTIKIIHD LSNIPQQDSI IDCFIYAFNP MDSFIPYYDI IDQKDCLTIK NSYELKVSES SKIERTKTL EEKSKPLVRP TARSKTTPEE TTGRKSKSKD MGLRSTALLA KMKKDRDDKE TSRQNELRKR KEENLQKINK QNPAREAQMK ELNNLFVEDD LDTEEVNGGS KPNSPKETDS NDKDKNNDL EDLLETTAED SLMDVPKIQQ TKPSETEHSK EPKSEEEPSS FIDEDGYIVT KRPATSTPPR KPSPVVKRAL SSSKKQETPS SNKRLKKQGT LESFFKRKAK
Source	Yeast
Target Names	POL32
Protein Names	Recommended name: DNA polymerase delta subunit 3
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.