



# Recombinant Mesocricetus auratus M-phase inducer phosphatase 3 (CDC25C)

<b>Product Code</b>	CSB-BP004996MRG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P48968
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mesocricetus auratus (Golden hamster)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	STGPFPSR REESSVSAPS FRFSQRKMLN LLLERNTSFT QDFPRSPGDK LLDSTNLSIL SGGTPKRCLD LSNLSNGEMS ASPLITSADF DDTGSLDSSG PQDVQLTEKN HHQDPMKGIP VQLLCSTPNA LDHSHRKKDA VRGLSANKEN INTNLKTLQW ESPRIPRFQN TPGDPLASPL PLLGNGVSM DTEVRS LGSPI TAVPKLSKNL NLEDQEEISE EPMEFSLEDH DTKECVLPTV SGKHQDLKYI TPDTVAALLS GKFQGLIEKF YIIDCRYPYE YLGGHILGAI NLCSQKELHE FFLKKPIVPL DIQKRVIIVF LCEFSSERGP RMCRLRRKD RALNQYPALY YPELYILKGG YRDFFPPEYTE LCEPQGYCPM HHQDHQAELL MWRNQSKAQE GERQLSEQIA LLMKKGVS LP
<b>Source</b>	Baculovirus
<b>Target Names</b>	CDC25C
<b>Protein Names</b>	Recommended name: M-phase inducer phosphatase 3 EC= 3.1.3.48 Alternative name(s): Dual specificity phosphatase Cdc25C
<b>Expression Region</b>	2-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.