



Recombinant 60S ribosomal protein L3 (rpl-3)

Product Code	CSB-MP020224CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50880
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSHRKFSAPR HGHMGFTPCK RSRTYRGRIK AFPKDDKSKP IHLTAFLGYK AGMTHIVRDV DKPGSKVNKK EVVEAVTIVE TPPMVIAGVT GYVDTPQGPR ALTTIWAHEL SEEARRRFYS NWAQSKKKAF TKYAKKWQDE DGKKLIEADF AKLKKYCSSI RVIAHTQMKI LRRRQKKAHL VEIQVNGGTI EQKVDWAREH LEKQVQVDTV FAQDEMIDI GVTRGHGFKG VTSRWHTKKL PRKTHKGLRK VACIGAWHPS RVAFTVARAG QKGFHHRTII NNKIYRIGKS ALTEEGKNNG STFDLTQKT ITPMGGFPRY GIVNQDYIML RGAVLGPKKR LITLRKSLIT QTKRVAHEKI NLKWIDTSSK TGHGRFQTTA EKRAFMGKLL RDFLAEEAEAK A
Source	Mammalian cell
Target Names	rpl-3
Protein Names	Recommended name: 60S ribosomal protein L3
Expression Region	1-401
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.