



Recombinant *Xenopus laevis* Retinoic acid receptor alpha (rara)

Product Code	CSB-EP019338XBE-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P51126
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MSSKDNTCPP PGPHINGFH VPHYAFFPH MLGGMSXTGG LPGVQHQPPL SGYSTPSPAT IETQSTSSEE IVSPPTPPP LPRIYKPCFV CQDKSSGYHY GVSACEGCKG FFRRSIQKNM VYTCHRDKNK IINKVTRNRC QYCRLQKCFE VGMSKESVRN DRNKKKKEP KPEAIESYIL SPETQDLIEK VQKAHQETFP ALCQLGKYTT SFSSEQRVSL DIDLWDFSE LSTKCIKTV EFAKQLPGFT TLTIADQITL LKSACLDILI LRICTRYTPD QDTMTFSDGL TLNRTQMHNA GFGPLTDLVF AFANQPVPLE MDDAETGLLS AICLICGDRQ DLEQPKVVDK LQEPLLEALK IYVRTRRPQK PHMFPKMLMK ITDLRTVSAK GAERVITLKM EIPGAMPLIQ EMLENSEGLD TLGGGASSDA PVTPVAPGSC SPSLSPSSTH SSPSTHSP
Source	E.coli
Target Names	rara
Protein Names	Recommended name: Retinoic acid receptor alpha Short name= RAR-alpha Alternative name(s): Nuclear receptor subfamily 1 group B member 1
Expression Region	1-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.