



Recombinant Danio rerio Fibroblast growth factor 3 (fgf3)

Product Code	CSB-EP008630DIL-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P48802
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	ES LAPRLTRTPR APCARGQACD PRQRRDAGGR GGVYEHLLGA PRRRKLYCAT KYHLQIHPNG KIDGSLEENN PLSILEITAV DVGVVAIKGL FSGRYLAMNE KGRLYASEVF NRECEFLERI HELGYNTYAS RHHATTQPPP TGSGIGGSKR RASSKRQWYV SINGKGRPRR GFKTRSTDKA SLFLPRVLAN KDHEMVRKLR ESQRHHTGSH RAPVGRAERR RRRHRGSKGH NRRADI
Source	E.coli
Target Names	fgf3
Protein Names	Recommended name: Fibroblast growth factor 3 Short name= FGF-3 Alternative name(s): Heparin-binding growth factor 3 Short name= HBGF-3
Expression Region	19-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.