



Recombinant Tubulin beta chain (TUBB1)

Product Code	CSB-EP025318EXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50259
Product Type	Recombinant Protein
Immunogen Species	Porphyra purpurea (Red seaweed) (Ulva purpurea)
Purity	>85% (SDS-PAGE)
Sequence	MREIVHIQAG QCGNQIGAKF WEVISEEHGI DSSGAYNGTS DLQLDRAEVY YNEGSGGRYV PRAVLVDLEP GVLDTIKAGP HGGLYRPDNF VAGQSGAGNN WAKGHYTEGA ELVDSVLDVV RREAEGCDCL QGFQVTHSLG GGTGSGMGTL LVSKIREEF DRMMCTYSVM PSPKVSDTVV EPYNCTLSVH QLVENADAVF CIDNEALYNI CYNTLKKPEP AYPDLNKLVS GVMGITSSL RFPGQLNSDL RKLAVNLVPF PRLHFFMIGY APLSADGVTA YRAKTVADLT KQMFDPKNMM ADCDPRNGRY LTASAYFRGH VSTKEVEEQM DAIQTKNSGQ FIDWIPNAIK ASVCDVAP TG ETMSAAFIGN STAIQDIFKR VGSHFSAMFK RKAFLHWYTG EGMDEMEFTE AESNMNDLVS ELQQYEAATV EGEEDDAYA EGAVVNGDQS YEDQYAA
Source	E.coli
Target Names	TUBB1
Protein Names	Recommended name: Tubulin beta chain Alternative name(s): Beta-tubulin
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.