



# Recombinant Rat N-acylglucosamine 2-epimerase (Renbp)

<b>Product Code</b>	CSB-MP019562RA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P51607
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMDPGLLVLQ DMEKERETLQ VWKQRVGQEL DSVIAFWMEH SHDQEHGGFF TCLGRDGQVY DHLKYVWLQG RQVWMYCRLY RTFERFRRVE LLDAAKAGGE FLLSYARVAP PGKKCAFVLT QDGRPVKVQR TIFSECFYTM AMNELWKVTG EMHYQREAVE MMDQIIHWVR EDPAGLGRPQ LSGTLATEPM AVPMMLLNLV EQLGEEDEEM TDKYAE LGDW CAHRILQH VQ RDGQVLE NV SEDGKELPGC LGRHQNPGHT LEAGWFLLQY ALRKGDPK LQ RHIIDK FLLL PFHSGWDPEH GGLFYFQDAD DLCPTQLEWN MKLWWPHT EA MIAFLMGYRD SGDPALLNLF YQVAEYTFHQ FRDPEYGEWF GYLNQEGKVA LTIKGGPFKG CFHVPRCLAM CEQILGALLQ RLG PAPLGS L PAVPTREGSK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Renbp
<b>Protein Names</b>	Recommended name: N-acylglucosamine 2-epimerase Short name= AGE EC= 5.1.3.8 Alternative name(s): GlcNAc 2-epimerase N-acetyl-D-glucosamine 2-epimerase Renin-binding protein Short name= RnBP
<b>Expression Region</b>	1-430
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	The gene product inhibits renin activity by forming a dimer with renin, a complex known as high molecular weight renin. The encoded protein contains a leucine zipper domain, which is essential for its dimerization with renin. The gene product can catalyze the interconversion of N-acetylglucosamine to N-acetylmannosamine, indicating that it is a GlcNAc 2-epimerase. Transcript variants utilizing alternative promoters have been described in the literature.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.