



Recombinant *Saccharomyces cerevisiae* Glycogenin-2 (GLG2)

Product Code	CSB-EP010071SVG-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P47011
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAKKVAICTL LYSRDYLPGA LTLAYQLQKL LKHAVVEDEI TLCLLIEKKL FGDEFKPEI ALIRSLFKEI IIIEPLKDQE KSIEKNKANL ELLKRPELSH TLLKARLWEL VQFDQVFLD ADTLPLNKEF FEILRLYPEQ TRFQIAAVPD IGWPDMFNTG VLLLIPDLDL ATSLQDFLIK TVSIDGADQG IFNQFFNPIC NYSKEVLHKV SPLMEWIRLP FTYNVTMPNY GYQSSPAMNF FQQHIRLIHF IGTFKPWSRN TTDYDDHYYQ LWRSTQRELY SECHLSNYFT HLQLGNIETE TNFYHEPPCL QDLLNHGKRE NQKHVDLDIT SVDRNASQKS TAEKHDIKPK TSKPQSAFKF DWESTDYLDLDR VQRAFPKPDT
Source	E.coli
Target Names	GLG2
Protein Names	Recommended name: Glycogenin-2 EC= 2.4.1.186 Alternative name(s): Glycogen synthesis initiator protein 2 Glycogenin glucosyltransferase 2
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.