



Recombinant *Saccharomyces cerevisiae* UBX domain-containing protein 6 (UBX6)

Product Code	CSB-MP025537SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P47049
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MYEMSGIDSL FHDRVVDHYS HTSEQVIVVY ISSAAGDNSW LHQWFKPGNL SDEERENILW VRLVNGTKEC LLFKSIFPSS SAPSINILQN GLLCSIQGN SLSREQDPWE TFINGLQSVF KGQVTKRKLKFSKSNEEYQRV KRMIQNDKLE RKYVFQNTND PQRKPQKWKQ LTVTDNVSYK SQKGFLAQNY CTLQLKLPNG YTISNTFFPPQ TKLHKVRMWL DYNKYDDGTP YLFHRNIPRV TLTRNDELKS LQELDLLPRS TLILEPLEAN NKTFDYMEQS SLLHKVYSGL TSFWAKEPEV DASSRRLGYQ RLGTVNSNSA NYSLQKLSSL DMVSDGGGGG GGDSMTPSAY TTPRMYPNSG TSQLRQNVSE LNLSSNNSAS NTKVRTLGYS NNGNN
Source	Mammalian cell
Target Names	UBX6
Protein Names	Recommended name: UBX domain-containing protein 6
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.