



Recombinant *Saccharomyces cerevisiae* Actin-related protein 3 (ARP3)

Product Code	CSB-MP001249SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P47117
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSYLNPAVV MDNGTGLTKL GFAGNDSPSW VFPTAIATAA PSNTKKSSGV GAPSAVSNEA SYFGNSTSAT NFNGATGGLL SNNLSGKRGT EDLDFYIGNE ALVASQGPSY SLSYPIRHGQ VENWDHMERF WENSIFKYLR TEPEDHFFLL TEPPLNPPEN REQVAEIFFE SFNCAGLYIA VQAVLALAAS WTSSKVTDERS LTGTVIDSGD GVTHVIPVAE GYVIGSAIKN IPIAGRITL FIQSLRERG EADTSLRTAE KIKQEYCYVC PDIVKEFNKF DKDPSKFAQF VVENQEKTRR KVVDIGYERF LAPEIFFNPE IASSDFLTPL PTVVDQTIQA CPIDVRKGLY NNIVLSGGST MFKDFGRRQL RDLKSIVNNR IAQSELLSGT KSTGVDVSVI SHRQQRNAVW FGGSLLAQTA EFKGYCHTKK DYEEYGPEIV RNFSLFNMV
Source	Mammalian cell
Target Names	ARP3
Protein Names	Recommended name: Actin-related protein 3 Alternative name(s): Actin-like protein ARP3 Short name= Actin-like protein 3
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.