



Recombinant *Drosophila melanogaster* Transcription initiation factor TFIID subunit 12 (Taf12)

Product Code	CSB-MP023079DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49905
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSDLFTTFDS NGVAEHLHH NHNSTSSASG LLHDPPMASP SQHSPMTNNS NSSSQNGGPV SGLGTGTGPI SGGSKSSNHT SSAAGSENTP MLTKPRLTEL VREVDTTTQL DEDVEELLQ IIDDFVEDTV KSTSAFAKHR KSNKIEVRDV QLHFERKYNM WIPGFGTDEL RPYKRAAVTE AHKQRLALIR KTIKKY
Source	Mammalian cell
Target Names	Taf12
Protein Names	Recommended name: Transcription initiation factor TFIID subunit 12 Alternative name(s): TAFII30 alpha Transcription initiation factor TFIID 28-alpha kDa/22 kDa subunits p28-alpha/p22
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.