



Recombinant *Saccharomyces cerevisiae* Ubiquitin carboxyl-terminal hydrolase 8 (UBP8)

Product Code	CSB-EP025752SVG-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50102
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSICPHIQQV FQNEKSKDGV LKTCNAARYI LNHSVPKEKF LNTMKCGTCH EINSGATFMC LQCGFCGCWN HSHFLSHSKQ IGHIFGINSN NGLLFCFKCE DYIGNIDLIN DAILAKYWDD VCTKTMVPSM ERRDGLSGLI NMGSTCFMSS ILQCLIHNPY FIRHSMSQIH SNNCKVRSPD KCFSCALDKI VHELYGALNT KQASSSSTST NRQTGFYLL TCAWKINQNL AGYSQQDAHE FWQFIINQIH QSYVLDLPNA KEVSRANNKQ CECIVHTVFE GSLESSIVCP GCQNNSKTTI DPFLDLSLDI KDKKKLYECL DSFHKKEQLK DFNYHCGECN STQDAIKQLG IHKLPVSLVL QLKRFHELLN GSNRKLDDFI EFPTYLNMKN YCSTKEKDKH SENGKVPDII YELIGIVSHK GTVNEGHYIA FCKISGGQWF KFNDSMVSSI SQEEVLKEQA YLLFYTIRQV N
Source	E.coli
Target Names	UBP8
Protein Names	Recommended name: Ubiquitin carboxyl-terminal hydrolase 8 EC= 3.4.19.12 Alternative name(s): Deubiquitinating enzyme 8 Ubiquitin thioesterase 8 Ubiquitin-specific-processing protease 8
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.