



Recombinant Oomycete-like sp. Tubulin beta-3 chain (TUBB3)

Product Code	CSB-BP025323OBR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50261
Product Type	Recombinant Protein
Immunogen Species	Oomycete-like sp. (strain MacKay2000)
Purity	>85% (SDS-PAGE)
Sequence	MREIVHIQAG QCGNQIGAKF WQVISDEHGV DPTGTYKGDS DLQLERINVY YNEATGGRYV PRAVLIDLEP GTMDAVRAGP FGQLFRPDNF VFGQTGAGNN WAKGHYTEGA ELIDSVLDVV RNEAESCDCD QGFQMTSLG GGTGSGMGTL LIAKLREEYP DRVMMTFSVC PSPKVSOTVV EPYNATLSVH QIVENTDES FVLDNEALYDI CFRTLKLTTP TYGDLNHLVC VCMSGVTSSL RFPGQLNCDL RKVAVNLIPF PRLHFFMVGF APLTSRGSQQ YRALTVPILT QQCFDAKNMM CAADPRHGRF LTASCMFRGR MSTKEVDEQM LNVQTKNSSY FVEWIPNNIK ASVCDIPPKG LKMSSTFVGN STAIQEMFKR VGEQFTAMYK RKAFLHWYTG EGMDDMEFTE AESNMNDLVS EYQQYQEASA DDEADEFDEE EGDEE
Source	Baculovirus
Target Names	TUBB3
Protein Names	Recommended name: Tubulin beta-3 chain Alternative name(s): Beta-3-tubulin
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.