



Recombinant Methylophilus methylotrophus Electron transfer flavoprotein subunit alpha (etfA)

Product Code	CSB-MP007844MFP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53571
Product Type	Recombinant Protein
Immunogen Species	Methylophilus methylotrophus (Bacterium W3A1)
Purity	>85% (SDS-PAGE)
Sequence	SKILVIAEH RRNDLRPVSL ELIGAANGLK KSGEDKVVVA VIGSQADAFV PALSVNGVDE LVVVKGSSID FDPDVFEASV SALIAAHNPS VVLLPHSVDS LGYASSLASK TGYGFATDVY IVEYQGDELV ATRGGYNQKV NVEVDFFPGKS TVVLTIRPSV FKPLEGAGSP VVSNVDAPSV QRSRQNKDYV EVGGGNDIDI TTVDFIMSIG RGIGEETNVE QFRELADEAG ATLCCSRPIA DAGWLPKSRQ VGQSGKVVGS CKLYVAMGIS GSIQHMAGMK HVPTIIAVNT DPGASIFTIA KYGIVADIFD IEEELKAQLA A
Source	Mammalian cell
Target Names	etfA
Protein Names	Recommended name: Electron transfer flavoprotein subunit alpha Short name= Alpha-ETF Alternative name(s): Electron transfer flavoprotein large subunit Short name= ETFLS
Expression Region	2-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.