



Recombinant Rat Vitamin K-dependent protein S (Pros1)

Product Code	CSB-EP018754RA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53813
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	ANTLLEETK KGNLERECIE ELCNKEEARE VFENNPETDY FYPKYLGCLG AFRVGAFSAA RQSANAYPDL RSCVNAIPDQ CDPMPCNEDG YLSCKDGGGA FTCICKPGWQ GDKCQFDINE CKDPSNINGG CSQTCDNTPG SYHCSCKIGF AMLTNKKDCK DVDECSLKPS VCGTAVCKNI PGDFECECPN GYRYDPSSKS CKDVDECSSEN TCAQLCVNYP GGYSCYCDGK KGFKLAQDQR SCEGIPVCLS LDLDKNYELL YLAEQFAGVV LYLKFRLPDI TRFSAEFDNR TYDSEGIILY AESLDHSNWL LIALREGKIE VQFKNEFSTQ ITTGGNVINN GIWNMVSVEE LDDSVSIKIA KEAVMNINKL GSLFKPTDGF LDTKIYFAGL PRKVESALIK PINPRLDGC I RGWNLMKQGA LGAKEIVEGK QNKHCFLTVE KGSYYPGSGI AQFSIDYNNV TNAEGWQINV TLNIRPSTGT GVMLALVSGD TVPFALSLVD SGSGTSQDIL VVENSVAAH LEAITLCSEQ PSQLKCNINR NGLELWTPVR KDVIYSKDLQ RQLAILDKTM KGTVATYLG VPDISFSATP VNAFYSGCME VNINGVQLDL DEAIKHNDI RAHSCPSVRK IQKNF
Source	E.coli
Target Names	Pros1
Protein Names	Recommended name: Vitamin K-dependent protein S
Expression Region	42-675
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a vitamin K-dependent plasma protein that functions as a cofactor for the anticoagulant protease, activated protein C (APC) to inhibit blood coagulation. It is found in plasma in both a free, functionally active form and also in an inactive form complexed with C4b-binding protein. Mutations in this gene result in autosomal dominant hereditary thrombophilia. An inactive pseudogene of this locus is located at an adjacent region on chromosome 3.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.