



Recombinant *Xenopus laevis* Paired box protein Pax-6 (pax6)

Product Code	CSB-EP017492XBE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P55864
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	<pre>MQNSHSGVNV LGGVFNVRP LPDSTRQKIV ELAHSGARPC DISRILQVSN GCVSKILGRY YETGSIRPRA IGGSKPRVAT PEVVNKIAHY KRECPSIFAW EIRDRLLESG VCTNDNIPSV SSINRVLRL ASDKQQMGSE GMYDKLRMLN GQTATWGSRP GWYPGTSVPG QPAQEGCQPQ EGVGENTNSI SSNGEDSDEA QMRLQLKRKL QRNRTSFTQE QIEALEKEFE RTHYPDV FAR ERLA AKIDLP EARIQVWFSN RRAKWRREEK LRNQRQASN TPSHIPSSS FSASVYQPIP QPTTPVSSFT SGSMLGRTDT ALTNSYSALP PMPSTMGNN LPMQPPVPSQ TSSYSCMLPT SPSVNGRSYD TYTPPHMQTH MNSQPMGTSG TTSTGLISPG VSPVQVPGS EPDMSQYWPR LQ</pre>
Source	E.coli
Target Names	pax6
Protein Names	Recommended name: Paired box protein Pax-6
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.