



Recombinant Triticum aestivum Protein disulfide-isomerase (PDI)

Product Code	CSB-EP017342TQN-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P52589
Product Type	Recombinant Protein
Immunogen Species	Triticum aestivum (Wheat)
Purity	≥85% (SDS-PAGE)
Sequence	EEAAA AEEAAAPEA VLTLHADNFD DAIKHPFIL VEFYAPWCGH CKSLAPEYEK AAQLLSKHDP AIVLAKVDAN DEKNKPLAGK YEVQGFPTLK IFRSGGKNIQ EYKGPREAEG IVEYLKKQVG PASKEIKAPE DATYLEDGKI HIVGVFTEFS GTEFTNFLEL AEKLRSDYDF GHTVHANHLP RGDAVERPL VRLFKPFDEL VVDSKDFDVS ALEKFIDASS TPKVVTFDKN PDNHPYLLKY FQSNAPKAML FLNFSTGPFE SFKSAYYGAV EEFSGKDVKF LIGDIEASQG AFQYNGLKED QAPLILIQDS DSKKFLKEQV EAGQIVAWLK DYFDGKLTPE RKSEPIPEAN NEPVKVVVAD NIHDVVFKSA KNVLIEFYAP WCGHCKKLAP ILDEAAATLQ SEEDVVIKI DATANDVPGE FDVQGYPTLY FVTPSGKKVS YEGGRTADEI VDYIKKNKET AGQAAAAATE KAAEPAATEP LKDEL
Source	E.coli
Target Names	PDI
Protein Names	Recommended name: Protein disulfide-isomerase Short name= PDI EC= 5.3.4.1
Expression Region	26-515
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.