



# Recombinant Takifugu rubripes Proto-oncogene c-Fos (fos)

<b>Product Code</b>	CSB-MP008790TKY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P53450
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Takifugu rubripes (Japanese pufferfish) (Fugu rubripes)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMFTSFNAEC DSSSRCSASP VGDNLYYPSP AGSYSSMGSP QSQDFTDLTA SSASFIPTVT AISTPDLQW MVQPLISSVA PSHRAHPYSP SPSYKRTVMR SAASKAHGKR SRVEQTTPEE EEKRRIRER NKQAAAKCRN RRRELTDTLQ AETDQLEDEK SSLQNDIANL LKEKERLEFI LAAHQPICKI PSQMDTDFSV VSMSPVHACL STTVSTQLQT SIPEATTVTS SHSTFTSTSN SIFSGSSDSL LSTATVSNSV VKMTDL DSSV LEESLDLLAK TEAETARSVP DVNLSNSLFA AQDWEPLHAT ISSSDFEPLC TPVVTCTPAC TTLTSSFVFT FPEAETFPTC GVAHRRRSNS NDQSSDSLSS PTLAL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	fos
<b>Protein Names</b>	Recommended name: Proto-oncogene c-Fos Alternative name(s): Cellular oncogene fos
<b>Expression Region</b>	1-376
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.