



Recombinant Acetyl-CoA acetyltransferase (phaA)

Product Code	CSB-MP001134PBZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P54810
Product Type	Recombinant Protein
Immunogen Species	Paracoccus denitrificans
Purity	≥85% (SDS-PAGE)
Sequence	MTKAVIVSAA RTPVGSFLGS FANLPAHELG AIVLKAVVER AGIDPSEVSE TILGQVLTAQ QGQNPARQAH IKVGLPRESA AWWINQVCGS GLRTVALAAQ QVLLGDARIV VAGGQESMSL APHAAYIAPG QKMGDMKMLD TMIKDGLWDA FNDYHMGTTA ENVAGKWEIS RAEQDQFAVA SQNKAEAAQK AGKFADEIVP VTIKSRKGET VVDADEYIRH GATLEAMENV RPAFSKEGTV TAGNASGLND GAAAVLVMTE DEAARRGLTP LARIASYATA GVDPQIMGTG PIPASRKALE KAGWSVGDLD LVEANFAFAA QIAVNRDMG WDPSIVNVNG GAIAIGHPIG ASGCRILNTL LFEMQRRDAK KGLATLCIGG GMGVALCLER P
Source	Mammalian cell
Target Names	phaA
Protein Names	Recommended name: Acetyl-CoA acetyltransferase EC= 2.3.1.9 Alternative name(s): Acetoacetyl-CoA thiolase
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.