



# Recombinant *Saccharomyces cerevisiae* Zinc finger protein ZPR1 (ZPR1)

<b>Product Code</b>	CSB-EP026617SVG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P53303
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSEQKEDLFK PVGEAAAEVE DESIAEQNKA NDGVKLTGAQ DAMGHPVQEI ESLCMNCGKN GTTRLLLSI PYFREIIMS FDCPHCGFKN CEIQPASQIQ EKGSRYVLKV ECREFNQRV IKSETATCKF VELDIEIPAK RGQLTTVEGL LSEMIDLSQ DQEMRKSIDE ALYKKIDDFI QKVKSINCE PNTIPITFIL DDPAGNSWIE YKPGEPQHKW SHTQYVRTDE QNVQVGIITR DQLEQRRREQ LKQLANRERN PSESVKVGSA NPQFLSDATD IENFNNEVQT FRASCPSTQ ECETHMKPVN IPHFKEVIIM STVCDHCGYK SNEVKTGGAI PDKGRRITLY CDDAADLSRD ILKSETCSMV IPELHLDIQE GTLGGRFTTL EGLLRQVYEE LESRIFTQTS DSMDEATKAR WVEFFAKLKE AIAGKVKFTV IMEDPLAGSY IQNVYAPDPD PNMTIEDYER TKEQNEDLGL SDIKVE
<b>Source</b>	E.coli
<b>Target Names</b>	ZPR1
<b>Protein Names</b>	Recommended name: Zinc finger protein ZPR1
<b>Expression Region</b>	1-486
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.