



Recombinant *Saccharomyces cerevisiae* Ubiquinone biosynthesis monooxygenase COQ6 (COQ6)

Product Code	CSB-MP005805SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53318
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MFFSKVMLTR RILVRGLATA KSSAPKLTDV LIVGGGPAGL TLAASIKNSP QLKDLKTTLV DMVDLKD KLS DFYNSPDYF TNRIVSVTPR SIHFLENNAG ATLMHDRIQS YDGLYVTDGC SKATLDLARD SMLCMIEIIN IQASLYNRIS QYDSKKDSID IIDNTKVVNI KHSDPNDPLS WPLVTLNNGE VYKTRLLVGA DGFNSPTRRF SQIPSRGWMY NAYGVVASMK LEYPPFKLRG WQRFLPTGPI AHLPM PENNA TLVWSSSERL SRRLLSLPPE SFTALINAAF VLEDADMNYY YRTLEDGSMD TDKLIEDIKF RTEEIYATLK DESDIDEIYP PRVVSIIKT RARFPLKLTH ADRYCTDRVA LVGDAAHTTH PLAGQQLNMG QTDVHGLVYA LEKAMERGLD IGSSLSLEPF WAERYPSNNV LLGMADKLFK LYHTNFPPVV ALRTFGLNLT NKIGPVKNMI IDTLGGNEK
Source	Mammalian cell
Target Names	COQ6
Protein Names	Recommended name: Ubiquinone biosynthesis monooxygenase COQ6 EC=1.14.13.-
Expression Region	1-479
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.