



Recombinant Anopheles gambiae 40S ribosomal protein S3a (RpS3A)

Product Code	CSB-BP020444BZL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P52813
Product Type	Recombinant Protein
Immunogen Species	Anopheles gambiae (African malaria mosquito)
Purity	>85% (SDS-PAGE)
Sequence	AVGKNKGVSV KGGKKGSKKK VVDPFTRKDW YDVKAPNMFK NRQVGKTLVN RTQGTKIASD GLKGRVFEVS LADLQNEPDA ERSFRKFKLV AESVNGRDVL TNFHGMALTT DKLRSVMKKW QTLIECSVDV KTTDGFMLRV FCIGFTIKDS MSQRKTCYAQ HSQIKNIRAK MTAIKREIT STDLKGVVEK LLPDSIAKDI EKACQVVYPL HDVYIRKVKV LKKPRFDLSS LMELHGDGGG KAAEVSTGAA SGVVVDRPEG YEPPVQASV
Source	Baculovirus
Target Names	RpS3A
Protein Names	Recommended name: 40S ribosomal protein S3a Alternative name(s): C3 protein
Expression Region	2-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.