



# Recombinant *Saccharomyces cerevisiae* tRNA-specific adenosine deaminase 1 (TAD1)

<b>Product Code</b>	CSB-EP001327SVG-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P53065
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVSCQGTRPC IVNLLTMPSE DKLGEIISTR VINEYSKLKS ACRPIIRPSG IREWTILAGV AAINRDGGAN KIEILSIATG VKALPDSELQ RSEGKILHDC HAEILALRGA NTVLLNRIQN YNPSSGDKFI QHNDEIPARF NLKENWELAL YISRLPCGDA SMSFLNDNCK NDDFIKIEDS DEFQYVDRSV KTLIRGRLNF NRRNVVRTKP GRYDSNITLS KSCSDKLLMK QRSSVLNCLN YELFEKPVFL KYIVIPNLED ETKHHLEQSF HTRLPLNDNE IKFLNCLKPF YDDKLDEEDV PGLMCSVKLF MDDFSTEEAI LNGVRNGFYT KSSKPLRKHC QSQVSRFAQW ELFKIRPEY EGISYLEFKS RQKKRSQLII AIKNILSPDG WIPTRTDDVK
<b>Source</b>	E.coli
<b>Target Names</b>	TAD1
<b>Protein Names</b>	Recommended name: tRNA-specific adenosine deaminase 1 EC= 3.5.4.- Alternative name(s): tRNA-specific adenosine-37 deaminase
<b>Expression Region</b>	1-400
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.