



Recombinant *Saccharomyces cerevisiae* Exosome complex component CSL4 (CSL4)

Product Code	CSB-MP007888SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53859
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MACNFQFPEI AYPGKLICPQ YGTENKDGED IIFNYVPGPG TKLIQYEHNG RTLEAITATL VGTVRCEEEK KTDQEEEREG TDQSTEEEEKS VDASPNDVTR RTVKNILVSV LPGTEKGRKT NKYANNDFAN NLPKEGDIVL TRVTRLSLQR ANVEILAVED KPSPIDSGIG SNGSGIVAAG GGSGAATFSV SQASSDLGET FRGIIRSQDV RSTDRDRVKV IECFKPGDIV RAQVLSLGDG TNYLTTARN DLGVVFARAA NGAGGLMYAT DWQMMTSPVT GATEKRKCAK PF
Source	Mammalian cell
Target Names	CSL4
Protein Names	Recommended name: Exosome complex component CSL4 Alternative name(s): CEP1 synthetic lethal protein 4
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.