



# Recombinant Rat Serine/threonine-protein phosphatase 5 (Ppp5c)

<b>Product Code</b>	CSB-MP018582RA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P53042
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AMAEGERTE CAEPPRDEPP AEGTLKRAEE LKTQANDYFK AKDYENAIKF YSQAIELNPS NAIYYGNRSL AYLRTCEYGY ALGDATRAIE LDKKYIKGYY RRAASNALG KFRAALRDYE TVVKVKPNDK DAKMKYQECS KIVKQKA FER AIAGDEHRRS VVDSL DIESM TIEDEYSGPK LEDGKVTITF MKDLMQWYKD QKKLHRKCAY QILVQVKEVL CKLSTLVETT LKETEKITVC GDTHGQFYDL LNIFELNGLP SETNPYIFNG DFVDRGSFSV EVILTLFGFK LLYPDHFHLL RGNHETDNMN QIYGFEGEVK AKYTAQMYEL FSEVFEWLPL AQCINGKVL I MHGGLFSEDG VTLDDIRKIE RNRQPPDSGP MCDLLWSDPQ PQNGRSVSKR GVSCQFGPDV TKAFLEENQL DYIIRSHEVK AEGYEVAHGG RCVTVFSAPN YCDQMGNKAS YIHLQGS DLR PQFHQFTAVP HPNVKPMAYA NTLLQLGMM
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Ppp5c
<b>Protein Names</b>	Recommended name: Serine/threonine-protein phosphatase 5 Short name= PP5 EC= 3.1.3.16 Alternative name(s): Protein phosphatase T Short name= PPT
<b>Expression Region</b>	2-499
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.