



Recombinant *Saccharomyces cerevisiae* Actin-related protein 2/3 complex subunit 2 (ARC35)

Product Code	CSB-YP002127SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53731
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MLHLQPQNLL IQKTLNEAIE ALRKGSP LTM DRIVSDFDYT TYHISNTAED KSILLLSVKT KAWVSVSECQ LDGSLTLLKF LADHYSSLGG VTIPSEVEPG YDYTLQITLA ELVQESILQL SVLKTIILSF PFELAISKFI ELSQQQPAPV EAEITGGEVA ANGDNLTFTI KYRDEENIFI KPSNDRVTII FETIFQDETD KIFGKVFLQE FVDARKRNRQ IQSAPQVLYS HEPPELKLRL YQPPKVAEQS RRFITFVLFV RHFQTKELQF HSICQLTLFR NYFHYYHIKCS KAYMHSRMRF RVDSFIKVLN RAKVDEDDEN DELSAEGRQQ ARRTFTGRKI VY
Source	Yeast
Target Names	ARC35
Protein Names	Recommended name: Actin-related protein 2/3 complex subunit 2 Alternative name(s): Arp2/3 complex 34 kDa subunit Short name= p34-ARC
Expression Region	1-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.