



Recombinant *Saccharomyces cerevisiae* Nucleoporin SEH1 (SEH1)

Product Code	CSB-MP020972SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53011
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	<p>MQPFDSGHDD LVHDEVYDFY GRHVATCSSD QHIKVFKLDK DTSNWELSDS WRAHDSSIVA IDWASPEYGR IIASASYDKT VKLWEEDPDQ EECSGRRWNK LCTLNDSKGS LYSVKFAPAH LGLKLAACLG N DGILRLYDAL EPSDLRSWTL TSEMKVLSIP PANHLQSDFC LSWCPSRFSP EKLAVSALEQ AIIYQRGKDG KLHVAAKLPG HKSLIRSISW APSIGRWYQL IATGCKDGRI RIFKITEKLS PLASEESLTN SNMFDNSADV DMDAQGRSDS NTEEKAELQS NLQVELLSEH DDHNGEVWSV SWNLTGTILS SAGDDGKVRL WKATYSNEFK CMSVITAQQ</p>
Source	Mammalian cell
Target Names	SEH1
Protein Names	Recommended name: Nucleoporin SEH1 Alternative name(s): Nuclear pore protein SEH1 SEC13 homolog 1
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.