



Recombinant Replication factor C subunit 4 (rfc-4)

Product Code	CSB-MP019591CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53016
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MEEPMEVDNK RPKVLTWTEK YRPKTLDDIA YQDEVVTMLK GALQGRDLPH LLFYGPPGTG KTSAAALAFCR QLFPKNIFHD RVLDLNASDE RGIAVVRQKI QSFSKSSLGH SHREDVLKLK IILDEVDAM TREAQAAMRR VIEDFSKTRR FILICNYVSR LIPPVVSRCA KFRFKSLPAE IQVQRLRTIC DAEGTPMSDD ELKQVMEYSE GDLRRVCTL QSLAPILKSG DDNARNCYLR GSSDSSLISN VCKSILTADV PQIIALTKDI TKSCGTGVAFI RRCFQQLMDE DVINDENIGV MGKLVATCEK RILDGCDLEN NLLDFLLTLR ETIQ
Source	Mammalian cell
Target Names	rfc-4
Protein Names	Recommended name: Replication factor C subunit 4 Alternative name(s): Activator 1 subunit 4
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.