



Recombinant Serinus canaria Transcription factor AP-1 (JUN)

Product Code	CSB-MP011972SYL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P54864
Product Type	Recombinant Protein
Immunogen Species	Serinus canaria (Island canary) (Fringilla canaria)
Purity	>85% (SDS-PAGE)
Sequence	MSAKMEPTFY EDALSAGFAP PESGGYGYNN AKVLKQNMTL NLSDPSSNLK PHLRNKNADI LTSPDVGLLK LASPELERLI IQSSNGLITT TPTPTQFLCP KNVTDEQEGF AEGFVRALAE LHNQNTLPSV TSAAQPVTSG MAPVSSMAGS TSFNLSLHSE PPVYANLSNF NPNALSSAPN YNANSMGYAP QHHINPQMPV QHPRLQALKE EPQTVPEMPG ETPPLSPIDM ESQERIKAER KRMRNRIAS KCRKRKLERI ARLEEKVKTL KQNSELAST ANMLREQVAQ LKQKVMNHVN SGCQLMLTQQ LQTF
Source	Mammalian cell
Target Names	JUN
Protein Names	Recommended name: Transcription factor AP-1 Alternative name(s): Proto-oncogene c-Jun
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.