



Recombinant *Macaca fascicularis* Cellular tumor antigen p53 (TP53)

Product Code	CSB-EP024077MOV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P56423
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Purity	>85% (SDS-PAGE)
Sequence	MEEPQSDPSI EPPLSQETFS DLWKLLPENN VLSPLPSQAV DDLMLSPDDL AQWLTEDPGP DEAPRMSEAA PPMAPTPAAP TPAAPAPAPS WPLSSSVPSQ KTYHGSYGFR LGFLHSGTAK SVTCTYSPDL NKMFCQLAKT CPVQLWVDST PPGSRVRAM AIYKQSQHMT EVVRRCPHHE RCSDSGLAP PQHLIRVEGN LRVEYSDDRN TFRHSVVVPY EPPEVGS DCT TIHYNMCNS SCMGGMNRRP ILTITLED SGNLLGRNSF EVRVCAC PGR DRRTEENFR KKGEPCHQLP PGSTKRALPN NTSSSPQPKK KPLDGEYFTL QIRGRERFEM FRELNEALEL KDAQAGKEPA GSRAHSSHLK SKKGQSTSRH KKFMFKTEGP DSD
Source	<i>E.coli</i>
Target Names	TP53
Protein Names	Recommended name: Cellular tumor antigen p53 Alternative name(s): Tumor suppressor p53
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.