



Recombinant Mouse Protein Wnt-11 (Wnt11)

Product Code	CSB-BP026131MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P48615
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	IKWLAL SKTPAALALN QTQHCKQLEG LVSAQVQLCR SNLELMRTIV HAARGAMKAC RRAFADMRWN CSSIELAPNY LLDLERTRE SAFVYALSAA TISHTIARAC TSGDLPGCSC GPVPGPEPPGP GNRWGGCADN LSYGLLMGAK FSDAPMKVKK TGSQANKLMR LHNSEVGRQA LRASLETCKK CHGVSGSCSI RTCWKGLQEL QDVAADLKTR YLSATKVVHR PMGTRKHLVP KDLDIRPVKD SELVYLQSSP DFCMKNEKVG SHGTQDRQCN KTSNGSDSCD LMCCGRGYNP YTDRVVERCH CKYHWCCYVT CRRCERTVER YVCK
Source	Baculovirus
Target Names	Wnt11
Protein Names	Recommended name: Protein Wnt-11
Expression Region	25-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It encodes a protein which shows 97%, 85%, and 63% amino acid identity with mouse, chicken, and Xenopus Wnt11 protein, respectively. This gene may play roles in the development of skeleton, kidney and lung, and is considered to be a plausible candidate gene for High Bone Mass Syndrome.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.