



Recombinant Pan troglodytes Growth hormone variant (GH2)

Product Code	CSB-MP009408EQV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P58757
Product Type	Recombinant Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Purity	>85% (SDS-PAGE)
Sequence	FPTI PLSRLFDNAM LRAHRLYQLA YDTYQEFEEA YILKEQKYSF LQNPQTSLCF SESIPTPSNR VKTQQKSNLE LLRISLLLIQ SWLEPVQLLR SVFANSLVYG ASDSNVYRHL KDLEEGIQTL MWRLEDGSPR TGQIFNQSYS KFDTKSHNDD ALLKNYGLLY CFRKDMDKVE TFLRIVQCRS VEGSCGF
Source	Mammalian cell
Target Names	GH2
Protein Names	Recommended name: Growth hormone variant Short name= GH-V Alternative name(s): Growth hormone 2 Placenta-specific growth hormone
Expression Region	27-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.