



Recombinant Rat cAMP-dependent protein kinase type I-beta regulatory subunit (Prkar1b)

Product Code	CSB-BP018696RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P81377
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	ASPSCFHSE DEDSLKGCEM YVQKHGIQQV LKECIVHLCV AKPDRPLRFL REHFEEKLEKE ENRQILARQK SNSQCDSHDE EISPTPPNPV VKARRRRGGV SAEVYTEEDA VSYVRKVIPK DYKTMTALAK AISKNVLFSSH LDDNERSDIF DAMFPVTHID GETVIQQGNE GDNFYVIDQG EVDVYVNGEW VTNISEGGSF GELALIYGTP RAATVKAKTD LKLWIGIDRDS YRRILMGSTL RKRKMYEEFL SKVSILESLE KWERLTVADA LEPVQFEDGE KIVVQGE PGD DFYIITEGTA SVLQRRSPNE EYVEVGR LGP SDYFGEIALL LNRPRAATVV ARGPLKCVKL DRPRFERVLG PCSEILKRNI QRYNSFISLT V
Source	Baculovirus
Target Names	Prkar1b
Protein Names	Recommended name: cAMP-dependent protein kinase type I-beta regulatory subunit
Expression Region	2-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.