



Recombinant *Drosophila melanogaster* Homeobox protein orthopedia (otp)

Product Code	CSB-MP017279DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P56672
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	<p>MLNNLGANVL GPKDCDLPKA ALTALHAGVN SFGQLPVGPN LADLGGGAGG GGSSGGVPVG GGVARLHISG GLCDNSNALN GGNGSSGNGN GNNNNGNGNN NNSMQQQDQK RHRTRFTPAQ LNELERCFSK THYPDIFMRE EIAMRIGLTE SRVQVWFQNR RAKWKKRKKK TNVFRTPGAL LPSHGLPPFG ANITNIAMGD GLCGTGMFVG DRWSVGVNPM TAGFGQLNQS SPLSSSLNSG LNSGINMGSA LGAGSYQHYP LNALGDSMMY QHSVGGVSCG PSGSPSATT PNMNSCSSVT PPPLSAQPNS SQNELNGEPM PLHQQQQQQT HQHQQQQTHQ HHPMAPPTPT QQQQQLPQSM QSPSDGANDT LHSIAALRRR ASELNAIPSY LQMAPHNYEH YNSNSNSVY</p>
Source	Mammalian cell
Target Names	otp
Protein Names	Recommended name: Homeobox protein orthopedia
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.