



Recombinant Saguinus oedipus Ras-related protein Ral-A (RALA)

Product Code	CSB-EP019296SWB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P63320
Product Type	Recombinant Protein
Immunogen Species	Saguinus oedipus (Cotton-top tamarin)
Purity	>85% (SDS-PAGE)
Sequence	MAANKPKGQN SLALHKVIMV GSGGVGKSAL TLQFMYDEFV EDYEPTKADS YRKKVLDGE EVQIDILDTA GQEDYAAIRD NYFRSGEGFL CVFSITEMES FAATADFREQ ILRVKEDENV PFLLVGNKSD LEDKRQVSVE EAKNRADQWN VNYVETSAKT RANVDKVFFD LMREIRARKM EDSKEKNGKK KRKSLAKRIR ERC
Source	E.coli
Target Names	RALA
Protein Names	Recommended name: Ras-related protein Ral-A
Expression Region	1-203
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.