



# Recombinant Mouse Inositol polyphosphate 1-phosphatase (Inpp1)

<b>Product Code</b>	CSB-MP011729MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P49442
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSDILLELLC VSEKAANIAR ACRQQETLFQ LLIEEKKGAE KKKFAADFK TLADVLVQEV IKQNMENKFP GLGKKVFGEE SNEFTNDLGE KITVELQSTE EETAELLSKV LNGNMPASEA LAQVVHEDVD LTDPTLESLD ISIPHESLGI WVDPIDSTYQ YIKGSANVKS NQGIFPSGLQ CVTILIGVYD LQTGLPLMGV INQPFASQNL TTLRWKGQCY WGLSYMGTNI HSLQLAISKS DSETQTENSD REFSSPFSAV ISTSEKDTIK AALSRVCGGS VFPAAGAGYK SLCVIQGLAD YIFSEDTTY KWDSCAAHAI LRAMGGGIVD MKECLERSPD TGLDLPQLLY HVENKGASGV ELWANKGGLI AYRSRNRLDT FLSRLIQNLG PVKTQA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Inpp1
<b>Protein Names</b>	Recommended name: Inositol polyphosphate 1-phosphatase Short name= IPP Short name= IPPase EC= 3.1.3.57
<b>Expression Region</b>	1-396
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene encodes the enzyme inositol polyphosphate-1-phosphatase, one of the enzymes involved in phosphatidylinositol signaling pathways. This enzyme removes the phosphate group at position 1 of the inositol ring from the polyphosphates inositol 1,4-bisphosphate and inositol 1,3,4-trisphosphate.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.