



# Recombinant Pan troglodytes General transcription factor IIH subunit 4 (GTF2H4)

<b>Product Code</b>	CSB-EP010013EQV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P60027
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pan troglodytes (Chimpanzee)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MESTPSRGLN RVHLQCRNLQ EFLGGLSPGV LDRLYGHPAT CLAVFRELP LAKNWVMRML FLEQPLPQAA VALWVKKEFS KAQEESTGLL SGLRIWHTQL LPGGLQLLIL NPIFRQNLRI ALLGGGKAWS DDTSQLGPKD HARDVPSLDK YAEERWEVVL HFMVGSPSAA VSQDLAQLLS QAGLMKSTEP GEPPCITSAG FQFLLDTPA QLWYFMLQYL QTAQSRGMDL VEILSFLFQL SFSTLGKDYS VEGMSDLLN FLQHREFGL VFQRKRKSRR YYPTRLAINL SSGVSGAGGT VHQPGFIVVE TNYRLYAYTE SELQIALIAL FSEMLYRFPN MVVAQVTRES VQQAIASGIT AQQIIHFLRT RAHPVMLKQT PVLPPITDQ IRLWELERDR LRFTEGVLYN QFLSQVDFEL LLAHARELGV LVFENSAKRL MVVTPAGHSD VKRFWKRQKH SS
<b>Source</b>	E.coli
<b>Target Names</b>	GTF2H4
<b>Protein Names</b>	Recommended name: General transcription factor IIH subunit 4 Alternative name(s): Basic transcription factor 2 52 kDa subunit Short name= BTF2 p52 General transcription factor IIH polypeptide 4 TFIIH basal transcription factor compl
<b>Expression Region</b>	1-462
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.