



Recombinant Escherichia coli Lipoyl synthase (lipA)

Product Code	CSB-EP012927ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P60716
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSKPIVMERG VKYRDADKMA LIPVKNVATE REALLRKPEW MKIKLPADST RIQGIKAAMR KNGLHSVCEE ASCPNLAECF NHGTATFMIL GAICTRRCPF CDVAHGRPVA PDANEPVKLA QTIADMALRY VVITSVDRDD LRDGGAQHFA DCITAIRESK PQIKIETLVP DFRGRMDRAL DILTATPPDV FNHNLENVPR IYRQVRPGAD YNWSLKLLER FKEAHPEIPT KSGLMVGLGE TNEEIIIVMR DLRRHGVTML TLGQYLQPSR HHLPVQRYVS PDEFDEMKA E ALAMGFTHAA CGPFVRSYH ADLQAKGMEV K
Source	E.coli
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.