



Recombinant Human 40S ribosomal protein S26 (RPS26)

Product Code	CSB-YP020407HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P62854
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	TKKRRNNGR AKKGRGHVQP IRCTNCARCV PKDKAIKKFV IRNIVEAAAV RDISEASVFD AYLPLKLYVK LHYCVSCAIH SKVVRNRSRE ARKDRTPPPR FRPAGAAPRP PPKPM
Source	Yeast
Target Names	RPS26
Protein Names	Recommended name: 40S ribosomal protein S26
Expression Region	2-115
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.