



# Recombinant *Macaca fuscata fuscata* Cellular tumor antigen p53 (TP53)

<b>Product Code</b>	CSB-EP024077MMX
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P61260
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Macaca fuscata fuscata</i> (Japanese macaque)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEEPQSDPSI EPPLSQETFS DLWKLLPENN VLSPLPSQAV DDLMLSPDDL AQWLTEDPGP DEAPRMSEAA PPMAPTPAAP TPAAPAPAPS WPLSSSVPSQ KTYHGSYGFR LGFLHSGTAK SVTCTYSPDL NKMFCQLAKT CPVQLWVDST PPGSRVRAM AIYKQSQHMT EVVRRCPHHE RCSDSGLAP PQHLIRVEGN LRVEYSDDRN TFRHSVVPY EPPEVGS DCT TIHYNMCNS SCMGGMNRRP ILTITLED SGNLLGRNSF EVRVCACPGR DRRTEENFR KKGEPCHQLP PGSTKRALPN NTSSSPQPKK KPLDGEYFTL QIRGRERFEM FRELNEALEL KDAQAGKEPA GSRAHSSHLK SKKGQSTSRH KKF MFKTEGP DSD
<b>Source</b>	<i>E.coli</i>
<b>Target Names</b>	TP53
<b>Protein Names</b>	Recommended name: Cellular tumor antigen p53 Alternative name(s): Tumor suppressor p53
<b>Expression Region</b>	1-393
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.