



Recombinant Danio rerio Zinc finger protein ZFPM1 (zfpm1)

Product Code	CSB-MP026463DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P60622
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSEMVHSRLK QGQGPAQAQQ SFYPPGSPAS VHKGATCFEC DITFNNINNF YVHKRLYCSS RHQQGETGGL VKEGAVTAAA PPVSHAASPQ ARPVSRAASA SPSCPDPAAG GTASEPKVVE VKIEDPGLKD ATCSSSSEGE GPGGGQASEG SQSPSGSAED QDDDPTRTFC QACNIRFSRH DNYIVHKRFY CASRHDPTNQ RPHSGKAAFL PQPIRTRKRK KMYE
Source	Mammalian cell
Target Names	zfpm1
Protein Names	Recommended name: Zinc finger protein ZFPM1 Alternative name(s): Friend of GATA protein 1 Short name= FOG-1 Short name= Friend of GATA 1 Zinc finger protein multitype 1
Expression Region	1-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.