



Recombinant Bordetella pertussis Single-stranded DNA-binding protein (ssb)

Product Code	CSB-MP022706BUA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P66846
Product Type	Recombinant Protein
Immunogen Species	Bordetella pertussis (strain Tohama I / ATCC BAA-589 / NCTC 13251)
Purity	>85% (SDS-PAGE)
Sequence	MASVNVKILV GNLGRDPEVR YSPDGAAICN VSIATTSQWK DKASGERREE TEWHRVVMYN RLAEIAGEYL KKGRSVYIEG RLKTRKWQDK DTGADRYSTE IVADQMQLG GRDSGGDSGG GYGGGYDDAP RQQRAPAQRP AAPQRPAPQ AAPAANLADM DDDIPF
Source	Mammalian cell
Target Names	ssb
Protein Names	Recommended name: Single-stranded DNA-binding protein Short name= SSB Alternative name(s): Helix-destabilizing protein
Expression Region	1-166
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.