



# Recombinant *Saccharomyces cerevisiae* Cytochrome b-c1 complex subunit 7 (QCR7)

<b>Product Code</b>	CSB-EP025664SVG-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P00128
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPQSFTSIAR IGDYILKSPV LSKLCVPVAN QFINLAGYKK LGLKFDDLIA EENPIMQTAL RRLPEDESYA RAYRIIRAHQ TELTHHLLPR NEWIKAQEDV PYLLPYILEA EAAAKEKDEL DNIEVSK
<b>Source</b>	<i>E.coli</i>
<b>Target Names</b>	QCR7
<b>Protein Names</b>	Recommended name: Cytochrome b-c1 complex subunit 7 Alternative name(s): Complex III subunit 7 Complex III subunit VII Ubiquinol-cytochrome c reductase c reductase complex 14 kDa protein
<b>Expression Region</b>	1-127
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.