



Recombinant Human Arylsulfatase L (ARSL)

Product Code	CSB-EP002146HU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P51690
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	DISASRPNI LLLMADDLGI GDIGCYGNNT MRTPNIDRLA EDGVKLTQHI SAASLCTPSR AAFLTGRYPV RSGMVSSIGY RVLQWTGASG GLPTNETTFA KILKEKGYAT GLIGKWHLGL NCESASDHCH HPLHHGFDHF YGMPFSLMGD CARWELSEKR VNLEQKLNFL FQVLALVALT LVAGKLTHLI PVSWMPIVWS ALSAVLLLAS SYFVGALIVH ADCFLMRNHT ITEQPMCFQR TTPLILQEVA SFLKRKNKHGP FLLFVSFLHV HIPLITMENF LGKSLHGLYG DNVEEMDWMV GRILDTLDVE GLSNSTLIYF TSDHGGSELEN QLGNTQYGGW NGIYKGGKGM GGWEGGIRVP GIFRWPGVLP AGRVIGEPTS LMDVFPTVVR LAGGEVPQDR VIDGQDLLPL LLGTAQHSDH EFLMHYCERF LHAARWHQRD RGTMWKVFHFV TPVFQPEGAG ACYGRKVCPC FGKVVHHPD PLLFDLSRDP SETHILTPAS EPVFYQVMER VQQAVWEHQR TLSPVPLQLD RLGNIWRPWL QPCCGPFPLC WCLREDDPQ
Source	E.coli
Target Names	ARSL
Protein Names	Recommended name: Arylsulfatase E Short name= ASE EC= 3.1.6.-
Expression Region	32-589
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.