



Recombinant Papio anubis Alpha-1-antitrypsin (SERPINA1)

Product Code	CSB-MP021053ERI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P01010
Product Type	Recombinant Protein
Immunogen Species	Papio anubis (Olive baboon)
Purity	≥85% (SDS-PAGE)
Sequence	EDPQG DAAQKTDTPP HDQNHPTLNK ITPSLAEFAF SLYRQLAHQS NSTNIFFSPV SIATAFAMLS LGTKADTHSE ILEGLNFNLT EIPEAQVHEG FQELLRTLNK PDSQLQLTTG NGLFLNKSLK VVDKFLEDVK NLYHSEAFSV NFEDTEEAKK QINNYVEKGT QGKVVDLVKE LDRDVFALV NYIFFK GKWE RPFVEATEE EDFHVDQATT VKVPMRRLG MFNIYHCEKL SSWVLLMKYL GNATAIFFLP DEGKLQHLEN ELTHDIITKF LENENRRSAN LHLPKLAITG TYDLKTVLGH LGITKVFSNG ADLSGVTEDA PLKLSKAVHK AVLTIDEKGT EAAGAMFLEA IPMSIPPEVK FNKPFVFLMI EQNTKSPLFI GKVVNPTQK
Source	Mammalian cell
Target Names	SERPINA1
Protein Names	Recommended name: Alpha-1-antitrypsin Short name= AAT Alternative name(s): Alpha-1 protease inhibitor Alpha-1-antiproteinase Serpin A1
Expression Region	16-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.